

心肌标志物复合非定值质控品 Plus 赋值表

(批号: 水平 3 202409001)

项目	仪器型号	方法学	单位	水平 3	
				靶值	靶值范围
脑钠肽(BNP)	ABBOTT AEROSET/ARCHITECT SYSTEMS	CMIA	pg/mL	1584	1263-1905
	ABBOTT ALINITY SYSTEMS	CMIA	pg/mL	1549	1239-1859
	BECKMAN COULTER ACCESS/DXI SYSTEMS	Chemiluminescence	pg/mL	1786	1429-2143
	SIEMENS ATELICA SYSTEMS	Chemiluminescence	pg/mL	1338	1071-1606
	MINDRAY CLIA SYSTEM	Chemiluminescence	pg/mL	1198	958-1438
	AUTOLUMO A6 SYSTEM	Chemiluminescence	pg/mL	1042	834-1250
肌钙蛋白 I(cTnI)	ABBOTT AEROSET/ARCHITECT SYSTEMS	CMIA	ng/mL	20.9	16.7-25.1
	ORTHO-CLINICAL VITROS SYSTEMS	Chemiluminescence	ng/mL	37.0	29.6-44.4
	TOSOH AIA SYSTEMS	EIA	ng/mL	40.7	28.5-52.9
	BECKMAN COULTER ACCESS/DXI SYSTEMS	Chemiluminescence	ng/mL	8.61	6.65-10.6
	MINDRAY CLIA SYSTEM	Chemiluminescence	ng/mL	16.277	13.0-19.5
	AUTOLUMO A6 SYSTEM	Chemiluminescence	ng/mL	0.692	0.554-0.830
超敏肌钙蛋白 I(hs-cTnI)	ABBOTT AEROSET/ARCHITECT SYSTEMS	CMIA	pg/mL	15119	11612-18626
	ABBOTT ALINITY SYSTEMS	CMIA (04Z21)	pg/mL	>3600	**
	ABBOTT ALINITY SYSTEMS	CMIA (08P13)	pg/mL	15577	12461-18692
	BECKMAN COULTER ACCESS/DXI SYSTEMS	Chemiluminescence	pg/mL	>27000	**
	BIOMERIEUX VIDAS SYSTEMS	ELFA	pg/mL	26223	20600-31846
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	pg/mL	10891	8713-13069
	SIEMENS ATELICA SYSTEMS	Chemiluminescence	pg/mL	11223	8978-13468
	SIEMENS DIMENSION SYSTEMS	LOCI	pg/mL	11508	8747-14270
	SIEMENS DIMENSION VISTA	LOCI	pg/mL	8707	6860-10555
超敏 C 反应蛋白(hsCRP)	ABBOTT ALINITY SYSTEMS	Immunoturbidimetric	mg/dL	1.21	0.97-1.45
	BECKMAN COULTER AU SERIES	Immunoturbidimetric	mg/dL	1.26	1.01-1.51
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	Immunoturbidimetric	mg/dL	1.32	1.06-1.58
C 反应蛋白(CRP)	ABBOTT AEROSET/ARCHITECT SYSTEMS	Immunoturbidimetric	mg/dL	1.22	0.98-1.46
	BECKMAN COULTER AU SERIES	Immunoturbidimetric	mg/dL	1.15	0.92-1.38
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	Immunoturbidimetric	mg/dL	1.11	0.89-1.33
	MINDRAY BS SYSTEM	Immunoturbidimetric	mg/L	11.9	9.52-14.28
同型半胱氨酸(HCY) *	ABBOTT AEROSET/ARCHITECT SYSTEMS	CMIA	μmol/L	31.5	25.2-37.8
	ABBOTT ALINITY SYSTEMS	CMIA	μmol/L	28.9	23.1-34.7
	BECKMAN COULTER AU SERIES	Immunoturbidimetric	μmol/L	35.0	28.0-42.0
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	Enzymatic	μmol/L	48.9	39.1->50.0
	SIEMENS IMMULITE SYSTEMS	Chemiluminescence	μmol/L	34.2	27.3-41.0
	MINDRAY BS SYSTEM	Enzymatic Cycling Method	μmol/L	32.1	25.68-38.52
D-二聚体(D-dimer)	BIOMERIEUX VIDAS SYSTEMS	ELFA	ng FEU/mL	9572	7657->10000
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	Tina-quant	ng FEU/mL	3460	2768-4152
	Snibe MAGLUMI X SYSTEM	Chemiluminescence	μg(FEU)/mL	46.8	37.44-56.16
	SUCCEEDER SF SYSTEM	Immunoturbidimetric	μg/mL	2.54	2.03-3.05
	WERFEN ACLTOP 50 SYSTEM	Immunoturbidimetric	μg/mL	12.27	9.82-14.72
	SYSMEX CS SYSTEM	Immunoturbidimetric	μg/mL	10.69	8.55-12.83
	BECKMAN COULTER AU SERIES	Immunoturbidimetric	μg/mL	6.30	5.04-7.56
肌酸激酶同工酶(CKMB)- 活性	ABBOTT AEROSET/ARCHITECT SYSTEMS	Immunoturbidimetric	U/L	134	107-161
	ABBOTT ALINITY SYSTEMS	Immunoturbidimetric	U/L	134	98.4-170
	BECKMAN COULTER AU SERIES	Immunoturbidimetric	U/L	59.7	47.7-71.6
	SIEMENS ADVIA SYSTEMS	Immunoturbidimetric, Enzymatic	U/L	86.6	69.2-104
	SIEMENS DIMENSION SYSTEMS	Immunoturbidimetric, Enzymatic	U/L	91.7	73.4-110
	MINDRAY BS SYSTEM	Immunosuppressive	U/L	73.0	58.4-87.6
肌酸激酶同工酶(CKMB)- 质量	BECKMAN COULTER ACCESS/DXI SYSTEMS	Chemiluminescence	ng/mL	106	84.8-127
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	ng/mL	77.9	62.4-93.5
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA (STAT)	ng/mL	75.2	60.2-90.2
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	ng/mL	93.6	71.8-115

	SIEMENS ATELICA SYSTEMS	Chemiluminescence	ng/mL	93.8	75.0-113
	SIEMENS DIMENSION SYSTEMS	EIA	ng/mL	95.9	76.7-115
	TOSOH AIA SYSTEMS	EIA	ng/mL	111	77.4-144
	ORTHO-CLINICAL VITROS SYSTEMS	CLIA	ng/mL	69.3	55.4-83.1
	MINDRAY CLIA SYSTEM	Chemiluminescence	ng/mL	99.6	79.7-119.5
肌红蛋白(MYO)	ABBOTT ALINITY SYSTEMS	CMIA	ng/mL	493	394-591
	ABBOTT ALINITY SYSTEMS	Immunoturbidimetric	ng/mL	377	301-452
	BECKMAN COULTER ACCESS/DXI SYSTEMS	Chemiluminescence	ng/mL	449	359-538
	BECKMAN COULTER AU SERIES	Immunoturbidimetric	ng/mL	380	304-456
	BIOMERIEUX VIDAS SYSTEMS	ELFA	ng/mL	638	510-765
	ORTHO-CLINICAL VITROS SYSTEMS	Chemiluminescence	ng/mL	763	610-915
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	ng/mL	321	256-385
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA (STAT)	ng/mL	337	269-404
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	ng/mL	458	366-550
	SIEMENS DIMENSION SYSTEMS	EIA	ng/mL	621	497-745
	SIEMENS DIMENSION VISTA	LOCI	ng/mL	430	344-516
	TOSOH AIA SYSTEMS	EIA	ng/mL	200	140-260
	MINDRAY CLIA SYSTEM	Chemiluminescence	ng/mL	1620	1134-2106
	Snibe MAGLUMI X SYSTEM	Chemiluminescence	ng/mL	456	319-593
	AUTOLUMO A6 SYSTEM	Chemiluminescence	ng/mL	267	187-347
N 端 B 型脑钠肽前体(NT-proBNP)	ABBOTT ALINITY SYSTEMS	CMIA	pg/mL	20800	16640-24959
	BIOMERIEUX VIDAS SYSTEMS	ELFA	pg/mL	24245	19396->25000
	ORTHO-CLINICAL VITROS SYSTEMS	CLIA	pg/mL	21470	17176-25764
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	pg/mL	26207	20966-31449
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA (STAT)	pg/mL	27732	22186-33278
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	pg/mL	28284	22628-33941
	SIEMENS DIMENSION SYSTEMS	LOCI	pg/mL	19101	14326-23876
	SIEMENS DIMENSION VISTA	LOCI	pg/mL	18001	14401-21601
	SIEMENS IMMULITE SYSTEMS	Chemiluminescence	pg/mL	> 35000	**
	Snibe MAGLUMI X SYSTEM	Chemiluminescence	pg/mL	30047	24038-36056
	Snibe MAGLUMI X SYSTEM	Chemiluminescence	pg/mL	24396	19517-29275
	AUTOLUMO A6 SYSTEM	Chemiluminescence	pg/mL	29974	23979-35969
肌钙蛋白 T(cTnT)	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA (STAT TNT-G5ST)	ng/mL	7.89	6.31-9.47
高敏肌钙蛋白 T(hs-cTnT)	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	pg/mL	9898	7918-> 10000
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA (STAT)	pg/mL	7222	4233-> 10000
降钙素原(PCT)	ABBOTT AEROSET/ARCHITECT SYSTEMS	CMIA	ng/mL	12.5	10.0-15.1
	ABBOTT ALINITY SYSTEMS	CMIA	ng/mL	12.2	9.75-14.6
	BIOMERIEUX VIDAS SYSTEMS	ELFA	ng/mL	13.3	10.6-16.0
	ORTHO-CLINICAL VITROS SYSTEMS	Chemiluminescence	ng/mL	11.2	8.94-13.4
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	ng/mL	10.1	8.07-12.1
	SIEMENS ATELICA SYSTEMS	Chemiluminescence	ng/mL	12.5	10.0-15.0
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	ng/mL	9.71	7.77-11.7
	SIEMENS DIMENSION SYSTEMS	Chemiluminescence	ng/mL	11.9	9.49-14.2
	Snibe MAGLUMI X SYSTEM	Chemiluminescence	ng/mL	12.1	9.68-14.5
	MINDRAY CLIA SYSTEM	Chemiluminescence	ng/mL	9.05	7.24-10.9
	AUTOLUMO A6 SYSTEM	Chemiluminescence	ng/mL	11.6	9.28-13.9
人绒毛膜促性腺激素(hCG)	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	mIU/mL	110	87.6-131
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	mIU/mL	130	104-156
	SIEMENS ATELICA SYSTEMS	Chemiluminescence	mIU/mL	137	110-165
β-人绒毛膜促性腺激素(β-hCG)	ABBOTT ALINITY SYSTEMS	CMIA	mIU/mL	145	116-174
	BECKMAN COULTER ACCESS/DXI SYSTEMS	Chemiluminescence	mIU/mL	156	125-187
	ORTHO-CLINICAL VITROS SYSTEMS	CLIA	mIU/mL	193	154-232
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	mIU/mL	144	115-173
	SIEMENS DIMENSION VISTA	LOCI	mIU/mL	200	160-241
髓过氧化物酶(MPO)	BECKMAN COULTER AU SERIES	Immunoturbidimetric	ng/mL	141	113-169

肌酸激酶(CK)	ABBOTT AEROSET/ARCHITECT SYSTEMS	NAC	U/L	572	458-686
	ABBOTT ALINITY SYSTEMS	NAC	U/L	612	490-735
	BECKMAN COULTER AU SERIES	Enzymatic, IFCC	U/L	483	313-654
	ORTHO-CLINICAL VITROS SYSTEMS	NAC, Colorimetric	U/L	439	351-527
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	NAC, IFCC	U/L	600	480-720
	SIEMENS ADVIA SYSTEMS	Enzymatic, IFCC	U/L	617	494-741
	SIEMENS ATELICA SYSTEMS	Enzymatic, IFCC	U/L	627	502-752
	SIEMENS DIMENSION SYSTEMS	Enzymatic, IFCC	U/L	607	486-728
	SIEMENS DIMENSION VISTA	Enzymatic, IFCC	U/L	568	454-682
和肽素(copeptin)*	未赋值	**	**	**	**
洋地黄毒苷(Digitoxin)	ABBOTT ALINITY SYSTEMS	Immunoturbidimetric	ng/mL	> 80.0	**
	ABBOTT AEROSET/ARCHITECT SYSTEMS	CEDIA	ng/mL	> 50.0	**
	ABBOTT AEROSET/ARCHITECT SYSTEMS	Latex Inhibition	ng/mL	> 80.0	**
	BECKMAN COULTER AU SERIES	Microparticles	ng/mL	79.8	63.8->90.0
	ROCHE COBAS C/E/HITACHI/MODULAR SYSTEMS	ECLIA	ng/mL	> 80.0	**
	SIEMENS ADVIA SYSTEMS	Chemiluminescence	ng/mL	70.8	56.7->80.0
	SIEMENS DIMENSION SYSTEMS	Magnetic Particle EIA	ng/mL	56.8	45.4-68.2
	SIEMENS DIMENSION VISTA	Magnetic Particle EIA	ng/mL	59.4	47.5-71.3
	THERMO FISHER SCIENTIFIC INDIKO SYSTEMS	CEDIA	ng/mL	> 50	**

注：以上靶值及范围具有批特异性，仅供本批质控品参考。因所用技术、仪器、试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据，根据良好的实验室规范要求，实验室须遵循相关技术规范确立自己的靶值和可接受范围。