

凝血八项非定值质控品赋值表
(批号：水平 1、水平 2、水平 3 202405231)

项目	仪器型号	单位	水平 1		水平 2		水平 3	
			靶值	范围	靶值	范围	靶值	范围
活化部分凝血活酶时间 (APTT)	Sysmex CS-5100	Seconds	30.9	24.7~37.1	67.9	54.3~81.4	95.8	76.6~115.0
	Werfen ACLTOP750	Seconds	34.5	27.6~41.3	55.0	44.0~66.0	65.7	52.6~78.8
	Mindray CX-6200	Seconds	28.1	22.5~33.7	94.5	75.6~113.4	111.5	89.2~133.8
	Mindray ExC810	Seconds	31.1	24.9~37.3	94.1	75.3~113.0	151.9	121.5~182.3
	Stago STA-R Evoiution	Seconds	34.8	27.9~41.8	79.7	63.7~95.6	94.7	75.8~113.7
	Rayto RAC2890	Seconds	35.0	28.0~42.0	73.9	59.1~88.7	111.4	89.1~133.7
	Succeeder SF-8100	Seconds	34.7	27.7~41.6	67.6	54.1~81.1	79.8	63.9~95.8
凝血酶原时间 (PT)	Sysmex CS-5100	Seconds	11.7	9.4~14.1	25.7	20.6~30.9	37.1	29.7~44.5
	Werfen ACLTOP750	Seconds	12.0	9.6~14.4	32.4	25.9~38.9	46.3	37.0~55.5
	Mindray CX-6200	Seconds	11.3	9.0~13.5	25.4	20.3~30.5	31.5	25.2~37.8
	Mindray ExC810	Seconds	14.5	11.6~17.4	25.9	20.7~31.1	31.5	25.2~37.8
	Stago STA-R Evoiution	Seconds	13.4	10.7~16.0	28.9	23.1~34.7	39.9	31.9~47.8
	Rayto RAC2890	Seconds	13.2	10.5~15.8	25.8	20.6~30.9	33.9	27.1~40.6
	Succeeder SF-8100	Seconds	12.8	10.2~15.3	30.1	24.1~36.2	44.6	35.7~53.5
凝血酶原时间 (PT) (INR)	Sysmex CS-5100	INR	1.00	0.80~1.20	2.29	1.83~2.75	3.36	2.69~4.03
	Werfen ACLTOP750	INR	1.09	0.87~1.31	3.05	2.44~3.66	4.39	3.51~5.27
	Mindray CX-6200	INR	0.96	0.77~1.15	2.25	1.80~2.70	2.82	2.26~3.38
	Mindray ExC810	INR	1.20	0.96~1.44	2.09	1.67~2.51	2.53	2.02~3.04
	Rayto RAC2890	INR	1.05	0.84~1.26	2.12	1.7~2.54	2.82	2.26~3.38
	Stago STA-R Evoiution	INR	1.01	0.81~1.21	2.62	2.1~3.14	3.90	3.12~4.68
凝血酶原时间 (PT) (%)	Sysmex CS-5100	%	92.7	74.2~111.2	25.3	20.3~30.4	15.7	12.6~18.8
	Werfen ACLTOP750	%	88.0	70.4~105.6	24.5	19.6~29.4	16.0	12.8~19.2
	Mindray CX-6200	%	104.1	83.3~124.9	34.4	27.5~41.2	26.2	20.9~31.4
	Mindray ExC810	%	71.2	57.0~85.5	26.9	21.6~32.3	21.4	17.1~25.7
	Stago STA-R Evoiution	%	99.7	79.7~119.6	27.7	22.1~33.2	18.3	14.7~22.0
凝血酶时间 (TT)	Sysmex CS-5100	Seconds	21.1	16.8~25.3	18.7	15.0~22.5	16.3	13.1~19.6
	Werfen ACLTOP750	Seconds	16.4	13.1~19.7	25.4	20.3~30.5	28.6	22.8~34.3
	Mindray CX-6200	Seconds	18.2	14.6~21.9	24.2	19.3~29.0	31.8	25.4~38.2
	Mindray ExC810	Seconds	17.0	13.6~20.4	20.9	16.7~25.1	25.7	20.5~30.8
	Stago STA-R Evoiution	Seconds	20.4	16.3~24.5	17.8	14.2~21.4	16.0	12.8~19.2
	Rayto RAC2890	Seconds	20.4	16.3~24.4	24.1	19.3~28.9	28.7	22.9~34.4
	Succeeder SF-8100	Seconds	18.7	14.9~22.4	17.9	14.3~21.5	17.3	13.9~20.8
抗凝血酶 III (AT-III)	Sysmex CS-5100	%	93.9	75.1~112.7	43.5	34.8~52.2	25.7	20.6~30.9
	Werfen ACLTOP750	%	93.5	74.8~112.2	40.0	32.0~48.0	28.0	22.4~33.6
	Mindray CX-6200	%	95.8	76.7~115.0	43.3	34.6~51.9	29.1	23.3~34.9
	Mindray ExC810	%	96.7	77.4~116.0	43.9	35.1~52.6	26.0	20.8~31.2
	Stago STA-R Evoiution	%	108.3	86.7~130.0	57.3	45.9~68.8	38.3	30.7~46.0
	Rayto RAC2890	%	112.2	89.8~134.7	49.8	39.8~59.7	20.0	16.0~24.0
D-二聚体 (D-Dimer) (FEU)	Sysmex CS-5100	mg/L	1.53	1.22~1.84	3.72	2.98~4.46	4.95	3.96~5.94
	Mindray CX-6200	ng/mL	1476.3	1181.0~1771.5	3085.9	2468.7~3703.1	4417.0	3533.6~5300.3
	Mindray ExC810	μg/mL	1.09	0.87~1.31	3.51	2.81~4.21	5.66	4.53~6.79
	Stago STA-R Evoiution	mg/L	0.82	0.66~0.98	1.57	1.26~1.88	2.10	1.68~2.52
	Rayto RAC2890	mg/L	1.10	0.88~1.32	2.85	2.28~3.42	3.55	2.84~4.26
	Succeeder SF-8100	μg/mL	0.77	0.62~0.92	1.68	1.34~2.02	1.72	1.38~2.06

D-二聚体（D-Dimer） （DDU）	Werfen ACLTOP750	ng/mL	473.0	378.4~567.6	932.0	745.6~1118.4	1163.5	930.8~1396.2
纤维蛋白（原）降解产物（FDP）	Sysmex CS-5100	mg/L	4.11	3.29~4.93	9.31	7.45~11.17	12.41	9.93~14.89
	Mindray CX-6200	μg/mL	4.74	3.79~5.69	7.30	5.84~8.76	9.29	7.43~11.15
	Mindray ExC810	μg/mL	4.64	3.71~5.57	7.45	5.96~8.94	8.56	6.85~10.27
	Stago STA-R Evoiution	μg/mL	4.00	3.20~4.80	6.02	4.82~7.22	8.10	6.48~9.72
	Rayto RAC2890	mg/L	4.18	3.34~5.02	8.72	6.98~10.46	11.16	8.93~13.39
纤维蛋白原（FIB）	Sysmex CS-5100	g/L	2.75	2.20~3.30	1.29	1.03~1.55	1.11	0.89~1.33
	Werfen ACLTOP750	g/L	2.74	2.19~3.29	1.25	1.00~1.50	1.12	0.90~1.34
	Mindray CX-6200	mg/dL	289.5	231.6~347.4	126.2	101.0~151.5	125.8	100.7~151
	Mindray ExC810	mg/dL	257.7	206.1~309.2	151.3	121.1~181.6	136.3	109.1~163.6
	Stago STA-R Evoiution	g/L	2.83	2.26~3.40	1.42	1.14~1.70	1.33	1.10~1.60
	Rayto RAC2890	g/L	2.76	2.21~3.31	1.37	1.10~1.64	1.11	0.89~1.33
	Succeeder SF-8100	g/L	3.24	2.59~3.89	1.60	1.28~1.92	1.43	1.14~1.72
纤溶酶原（PLG）	Werfen ACLTOP750	%	87.5	70.0~105.0	37.5	30.0~45.0	24.0	19.2~28.8

注：以上靶值及范围具有批特异性，仅供本批质控品参考。因所用技术、仪器、试剂的不同，或因制造商检测方法的改变，均可导致实验室实际测得的数据偏离此报告所提供的数据，根据良好的实验室规范要求，实验室须遵循相关技术规范确立自己的靶值和可接受范围。